

¹H-NMR

¹H NMR (400 MHz, DMSO-d₆) δ 9.63 (s, 1H), 8.64 (s, 1H), 8.53 (dd, *J* = 7.9, 1.7 Hz, 1H), 7.69 (dd, *J* = 4.9, 1.7 Hz, 1H), 7.61 – 7.53 (m, 2H), 7.44 (dd, *J* = 7.8, 1.7 Hz, 1H), 7.34 – 7.28 (m, 2H), 7.24 (td, *J* = 7.5, 1.8 Hz, 1H), 7.17 (td, *J* = 7.6, 1.5 Hz, 1H), 7.07 (dd, *J* = 7.9, 4.9 Hz, 1H), 6.92 (dd, *J* = 7.9, 1.5 Hz, 1H), 1.33 (s, 9H).

