

¹H-NMR

¹H NMR (400 MHz, Chloroform-d) δ 9.39 – 9.24 (m, 1H), 9.02 (dd, *J* = 4.4, 1.8 Hz, 1H), 8.20 (dd, *J* = 8.0, 1.8 Hz, 1H), 7.92 (dd, *J* = 7.5, 1.6 Hz, 1H), 7.83 (d, *J* = 8.8 Hz, 1H), 7.79 – 7.66 (m, 3H), 7.54 (dd, *J* = 8.0, 4.4 Hz, 1H).

