

¹H-NMR

¹H NMR (400 MHz, Chloroform-d) δ 7.36 – 7.27 (m, 5H), 6.93 (d, *J* = 7.2 Hz, 1H), 6.80 – 6.74 (m, 2H), 6.71 (tt, *J* = 8.9, 2.3 Hz, 1H), 6.18 (d, *J* = 7.2 Hz, 1H), 5.36 (d, *J* = 7.2 Hz, 1H), 4.53 (p, *J* = 7.0 Hz, 1H), 3.47 (s, 2H), 1.38 (s, 12H).

