

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 4.96 (d, *J* = 2.7 Hz, 1H), 4.70 (s, 4H), 4.45 (t, *J* = 5.9 Hz, 1H), 4.35 (dd, *J* = 6.0, 2.6 Hz, 1H), 3.84 (d, *J* = 10.0 Hz, 1H), 3.76 (d, *J* = 3.2 Hz, 1H), 3.68 (d, *J* = 5.7 Hz, 1H), 3.61 – 3.44 (m, 5H), 3.44 – 3.13 (m, 13H), 3.09 (t, *J* = 9.3 Hz, 1H), 2.98 (t, *J* = 9.2 Hz, 1H), 2.86 (q, *J* = 5.5 Hz, 1H), 2.43 (dd, *J* = 11.8, 9.3, 4.0 Hz, 1H), 2.26 (s, 3H), 2.25 – 2.17 (m, 1H), 1.87 (dt, *J* = 12.9, 4.0 Hz, 1H), 0.84 (t, *J* = 12.1 Hz, 1H).

