

¹H-NMR

¹H NMR (400 MHz, Chloroform-*d*) δ 7.23 – 7.02 (m, 7H), 6.82 (t, *J* = 9.0 Hz, 2H), 6.75 (d, *J* = 8.8 Hz, 1H), 6.71 (d, *J* = 8.6 Hz, 1H), 6.47 (d, *J* = 8.7 Hz, 2H), 4.34 – 4.19 (m, 1H), 4.11 – 4.03 (m, 1H), 3.01 (d, *J* = 33.2 Hz, 2H), 2.71 – 2.39 (m, 8H), 0.92 (t, *J* = 7.4 Hz, 3H).

