

¹H-NMR

¹H NMR (400 MHz, Deuterium Oxide) δ 5.46 (d, *J* = 3.9 Hz, 1H), 5.01 (d, *J* = 3.8 Hz, 1H), 3.93 (ddd, *J* = 10.8, 8.2, 3.1 Hz, 1H), 3.89 – 3.83 (m, 1H), 3.79 – 3.58 (m, 5H), 3.55 (dd, *J* = 9.9, 3.9 Hz, 1H), 3.49 – 3.24 (m, 5H), 3.19 – 2.96 (m, 4H), 2.07 (dt, *J* = 12.9, 4.2 Hz, 1H), 1.38 (q, *J* = 12.5 Hz, 1H).

