

¹H-NMR

¹H NMR (400 MHz, Deuterium Oxide) δ 3.81 (ddd, *J* = 8.2, 6.0, 3.8 Hz, 1H), 3.74 (dd, *J* = 11.6, 2.9 Hz, 1H), 3.71 – 3.64 (m, 2H), 3.60 – 3.52 (m, 2H), 2.62 (dd, *J* = 12.7, 3.8 Hz, 1H), 2.55 (dd, *J* = 12.7, 8.2 Hz, 1H), 2.26 (s, 3H).

