

1H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 11.98 (s, 1H), 9.36 (s, 1H), 5.24 (s, 2H), 5.13 (tt, *J* = 5.5, 2.9 Hz, 1H), 3.69 (s, 3H), 3.29 (d, *J* = 6.9 Hz, 2H), 2.29 – 2.21 (m, 2H), 2.20 – 2.13 (m, 2H), 2.08 (s, 3H), 1.73 (d, *J* = 1.3 Hz, 3H).

