

¹H-NMR

¹H NMR (400 MHz, DMSO-d₆) δ 8.74 (dd, *J* = 4.2, 1.7 Hz, 1H), 8.26 (ddd, *J* = 8.3, 1.8, 0.8 Hz, 1H), 7.93 (d, *J* = 9.1 Hz, 1H), 7.48 (dd, *J* = 8.3, 4.2 Hz, 1H), 7.41 (dd, *J* = 9.1, 2.9 Hz, 1H), 7.37 (d, *J* = 2.8 Hz, 1H), 3.91 (s, 3H).

