

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 10.71 (s, 1H), 8.70 (s, 1H), 8.63 (s, 1H), 6.19 (d, *J* = 1.2 Hz, 1H), 5.72 (dd, *J* = 5.9, 1.2 Hz, 1H), 4.98 (ddd, *J* = 9.8, 5.9, 2.6 Hz, 1H), 4.07 (ddd, *J* = 20.2, 9.0, 4.8 Hz, 1H), 3.99 (ddd, *J* = 11.1, 9.0, 2.4 Hz, 1H), 3.88 (td, *J* = 10.0, 4.8 Hz, 1H), 2.55 (t, *J* = 7.3 Hz, 2H), 2.41 (td, *J* = 7.3, 1.7 Hz, 2H), 1.61 (dq, *J* = 18.0, 7.3 Hz, 4H), 0.93 (q, *J* = 7.5 Hz, 6H), 0.86 (d, *J* = 7.6 Hz, 1H).

