

¹H-NMR

¹H NMR (400 MHz, Chloroform-d) δ 8.73 (s, 1H), 8.56 (s, 1H), 7.68 (d, *J* = 37.9 Hz, 5H), 7.45 (s, 2H), 6.79 (d, *J* = 15.6 Hz, 1H), 6.62 (d, *J* = 10.1 Hz, 2H), 4.66 (s, 2H), 4.13 (d, *J* = 17.8 Hz, 3H), 3.87 (d, *J* = 1.8 Hz, 6H), 3.76 (s, 2H), 2.82 (s, 2H).

