


PROTON CDC13 E10026



Current Data Parameters
NAME E10026
EXPNO 10
PROCNO 1

F2 - Acquisition Parameters
Date_ 20191112
Time 11.03 h
INSTRUM spect
PROBHD Z116098_0169 (
PULPROG zg30
TD 65536
SOLVENT CDCl3
NS 16
DS 2
SWH 8012.820 Hz
FIDRES 0.244532 Hz
AQ 4.0894465 sec
RG 111.15
DW 62.400 usec
DE 6.50 usec
TE 299.2 K
D1 1.00000000 sec
TD0 1
SF01 400.1324708 MHz
NUC1 1H
P1 9.55 usec
PLW1 17.00000000 W

F2 - Processing parameters
SI 65536
SF 400.1300081 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00