



# Specification Sheet

Beijing psaitong Biotechnology Co., Ltd.  
Tel/Fax:010-60605840  
E-mail: [psaitong@jm-bio.com](mailto:psaitong@jm-bio.com)  
Web: [www.jm-bio.com](http://www.jm-bio.com)

---

Product Name: Methacrylated Alginate-甲基丙烯酯含量:0.20-0.30mmol/g  
Product Number: M14193  
CAS:  
Formula:  
Molecular Weight:  
Storage at: -20°C

Test	Specification
Appearance	White spongy
Methacryloyl ester	0.2-0.3mmol/g
pH(5% in water, 40° C)	6 - 8
Solubility in water(5%, 40° C, Shaker)	0-5Hour