

产品说明书

26 版

基本信息

产品编号:	N10034
产品名称:	3-氧代-N-[(3S)-2-氧代四氢呋喃-3-基]辛酰胺

产品简介:

一种促进转录激活因子表达的试剂。N-3-oxo-octanoyl-L-Homoserine lactone promotes the expression of the transcriptional activator (and LuxR homolog) TraR.

体外研究:

N-3-Oxo-octanoyl-L-homoserine lactone (3OC8-HSL) pretreatment (10 μ M; for 2 days) exerts strong inhibitory effect on pathogen proliferation in the leaves.

Plants treated with N-3-Oxo-octanoyl-L-homoserine lactone shows increased activities of defense-related enzymes including peroxidase, catalase, phenylalanine ammonialyase, and superoxide dismutase.

使用方法（供参考）:

一：体外实验

溶于乙醇（10mg/ml）、DMSO（30mg/ml）、DMF（30mg/ml）和氯仿（50mg/ml）。

浓度	质量 溶液体积	1mg	5mg	10mg
1 mM		4.1446 mL	20.7228 mL	41.4456 mL
5 mM		0.8289 mL	4.1446 mL	8.2891 mL
10 mM		0.4145 mL	2.0723 mL	4.1446 mL

储备液的保存方式：-80° C, 6 months; -20° C, 1 month.

二：体内实验

1、请依序添加每种溶剂：10% DMSO→40% PEG300→5% Tween-80→45% Saline

Solubility: ≥ 2.5 mg/mL (10.36 mM); 澄清溶液

此方案可获得 ≥ 2.5 mg/mL (饱和度未知) 的澄清溶液。

以 1 mL 工作液为例，取 100 μ L 25.0 mg/mL 的澄清 DMSO 储备液加到 400 μ L PEG300 中，混合均匀；再向上述体系中加入 50 μ L Tween-80，混合均匀；然后再继续加入 450 μ L 生理盐水 定容至 1 mL。

2、请依序添加每种溶剂：10% DMSO→90% (20% SBE- β -CD in Saline)

Solubility: ≥ 2.5 mg/mL (10.36 mM); 澄清溶液

此方案可获得 ≥ 2.5 mg/mL (饱和度未知) 的澄清溶液。

以 1 mL 工作液为例，取 100 μ L 25.0 mg/mL 的澄清 DMSO 储备液加到 900 μ L 20% 的 SBE- β -CD 生理盐水水溶液中，混合均匀。

2 g SBE- β -CD (磺丁基醚 β -环糊精) 粉末定容于 10 mL 的生理盐水中，完全溶解至澄清透明。

3、请依序添加每种溶剂：10% DMSO→90% Corn Oil

Solubility: ≥ 2.5 mg/mL (10.36 mM); 澄清溶液

此方案可获得 ≥ 2.5 mg/mL (饱和度未知) 的澄清溶液，此方案实验周期在半个月以上的动物实验酌情使用。

以 1 mL 工作液为例，取 100 μ L 25.0 mg/mL 的澄清 DMSO 储备液加到 900 μ L 玉米油中，混合均匀。



注：工作液建议您现用现配，当天使用。

<1mg/ml 表示微溶或不溶。

普西唐提供的所有化合物浓度为内部测试所得，实际溶液度可能与公布值有所偏差，属于正常的批间细微差异现象。

请根据产品在不同溶剂中的溶解度选择合适的溶剂配制储备液；一旦配成溶液，请分装保存，避免反复冻融造成的产品失效。

注意事项：

1、为了您的安全和健康，请穿实验服并戴一次性手套操作。

2、以上信息仅做参考交流之用。

储存： -20℃

