

产品说明书

26 版

基本信息

产品编号:	P10077
产品名称:	DL-苯丙氨酸

产品简介:

一种苯丙氨酸衍生物。

体外研究:

Amino acids and amino acid derivatives have been commercially used as ergogenic supplements. They influence the secretion of anabolic hormones, supply of fuel during exercise, mental performance during stress related tasks and prevent exercise induced muscle damage. They are recognized to be beneficial as ergogenic dietary substances.

使用方法（供参考）:

一：体外实验			
Water: 溶于水			
浓度 \ 质量 溶液体积	1mg	5mg	10mg
1 mM	6.0536 mL	30.2682 mL	60.5364 mL
5 mM	1.2107 mL	6.0536 mL	12.1073 mL
10 mM	0.6054 mL	3.0268 mL	6.0536 mL
储备液的保存方式: -80° C, 6 months; -20° C, 1 month.			

注意事项:

- 1、为了您的安全和健康，请穿实验服并戴一次性手套操作。
- 2、以上信息仅做参考交流之用。

储存: 室温

