



Specification Sheet

Product Name:	Methacrylated Alginate
Product Number:	M14193
CAS:	
Formula:	
Molecular Weight:	
Storage at:	-20°C

Test	Specification
Appearance	White spongy
Methacryloyl ester	0.2-0.3mmol/g
pH(5% in water,40°C)	6 - 8
Solubility in water(5%,40°C,Shaker)	0-5Hour

真实

准确

稳定

负责

了解更多
扫一扫



北京普西唐生物科技有限公司
天津市北辰区青光镇优谷新科园145栋
010-60605840 www.jm-bio.com
psaitong@jm-bio.com

*更多产品~请官网查询。
*本公司产品仅用于科研。