



Specification Sheet

Product Name:	Squalene
Product Number:	S80147
CAS:	7683-64-9
Formula:	C ₃₀ H ₅₀
Molecular Weight:	410.72
Storage at:	2-8℃

Test	Specification
Appearance	Colorless to light yellow clear Liquid
Acid value	≤1 mg KOH/g
Indine value	350 - 450
Peroxide value	≤7 meq/kg
Saponification value	≤7 mg KOH/g
Purity(HPLC)	≥98%
¹ H-NMR(CDCl ₃)	Conforms to Structure

