



## Specification Sheet

|                   |                                                               |
|-------------------|---------------------------------------------------------------|
| Product Name:     | U 93631                                                       |
| Product Number:   | U10010                                                        |
| CAS:              | 152273-12-6                                                   |
| Formula:          | C <sub>17</sub> H <sub>21</sub> N <sub>3</sub> O <sub>2</sub> |
| Molecular Weight: | 299.37                                                        |
| Storage at:       | -20℃                                                          |

| Test               | Specification            |
|--------------------|--------------------------|
| Appearance         | White to off-white solid |
| Purity(HPLC)       | ≥98%                     |
| <sup>1</sup> H-NMR | Conforms to Structure    |

