

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 9.05 (ddd, *J* = 24.2, 4.6, 2.1 Hz, 1H), 8.54 – 8.10 (m, 3H), 7.62 (ddd, *J* = 19.4, 8.1, 4.5 Hz, 2H), 7.44 (dd, *J* = 8.2, 2.4 Hz, 1H), 7.40 – 6.81 (m, 8H), 5.22 (dq, *J* = 57.5, 6.8 Hz, 1H), 4.99 – 4.81 (m, 1H), 4.66 – 4.42 (m, 1H), 4.17 – 3.98 (m, 2H), 3.80 – 3.48 (m, 1H), 2.28 (d, *J* = 16.6 Hz, 1H), 1.59 – 1.18 (m, 6H).

