

## <sup>1</sup>H-NMR

<sup>1</sup>H NMR (400 MHz, DMSO-d<sub>6</sub>) δ 8.09 (d, *J* = 7.8 Hz, 1H), 7.57 (d, *J* = 1.1 Hz, 1H), 6.80 (d, *J* = 1.1 Hz, 1H), 4.39 (td, *J* = 8.1, 5.2 Hz, 1H), 2.92 (ddd, *J* = 14.7, 5.3, 0.9 Hz, 1H), 2.81 (dd, *J* = 14.8, 8.4 Hz, 1H), 1.81 (s, 3H).

