

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 7.96 (d, *J* = 8.0 Hz, 1H), 7.38 – 7.17 (m, 5H), 3.96 (q, *J* = 7.1, 5.5 Hz, 1H), 3.84 (d, *J* = 3.0 Hz, 1H), 3.50 (ddd, *J* = 8.9, 6.1, 3.0 Hz, 1H), 2.91 (dd, *J* = 13.2, 8.7 Hz, 1H), 2.80 (dd, *J* = 13.1, 5.7 Hz, 1H), 1.68 – 1.49 (m, 3H), 0.85 (dd, *J* = 12.2, 5.6 Hz, 6H).

