

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 11.43 (s, 1H), 8.75 (d, *J* = 2.2 Hz, 1H), 8.37 (dd, *J* = 4.7, 1.5 Hz, 1H), 8.07 – 7.98 (m, 2H), 7.68 (s, 1H), 7.53 (d, *J* = 16.7 Hz, 1H), 7.36 (dd, *J* = 8.0, 4.7 Hz, 1H), 7.21 (dd, *J* = 9.9, 2.4 Hz, 1H), 7.10 (d, *J* = 16.7 Hz, 1H), 6.99 (ddd, *J* = 9.7, 8.8, 2.4 Hz, 1H).

