

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 7.29 – 7.16 (m, 3H), 7.05 (d, *J* = 8.2 Hz, 1H), 6.92 (t, *J* = 2.0 Hz, 1H), 6.79 (t, *J* = 8.6 Hz, 2H), 5.56 (s, 2H), 4.32 (tq, *J* = 6.7, 3.1 Hz, 4H), 3.27 (s, 4H), 3.19 (td, *J* = 6.2, 3.0 Hz, 4H), 1.78 (p, *J* = 5.8 Hz, 2H).

