

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 9.25 (br, 1H), 9.16 (br, 1H), 7.10 (d, *J* = 8.5 Hz, 2H), 6.86 (d, *J* = 8.2 Hz, 1H), 6.71 (d, *J* = 8.5 Hz, 2H), 6.28 (dd, *J* = 8.2, 2.4 Hz, 1H), 6.18 (d, *J* = 2.4 Hz, 1H), 4.14 (ddd, *J* = 10.5, 3.6, 1.9 Hz, 1H), 3.89 (t, *J* = 10.4 Hz, 1H), 3.01 (td, *J* = 10.1, 5.6, 3.6 Hz, 1H), 2.89 – 2.71 (m, 2H).

