

1H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 8.63 (dd, *J* = 5.2, 0.8 Hz, 1H), 8.26 (s, 1H), 7.97 (d, *J* = 8.2 Hz, 2H), 7.85 – 7.81 (m, 1H), 7.77 (dd, *J* = 5.1, 1.8 Hz, 1H), 7.30 (d, *J* = 7.9 Hz, 2H), 2.82 (dd, *J* = 8.4, 6.8 Hz, 2H), 2.36 (s, 3H), 1.82 – 1.70 (m, 2H), 0.94 (t, *J* = 7.3 Hz, 3H).

