

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 8.14 (d, *J* = 2.2 Hz, 1H), 7.61 (d, *J* = 2.2 Hz, 1H), 7.58 (dd, *J* = 6.1, 3.1 Hz, 1H), 7.32 (d, *J* = 3.0 Hz, 1H), 7.31 (s, 1H), 5.69 (s, 2H), 4.50 (s, 2H), 4.33 (s, 3H), 4.05 (s, 7H), 2.61 (s, 3H).

