

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 9.43 (s, 1H), 9.16 (s, 1H), 8.45 (s, 1H), 7.82 (s, 1H), 7.79 (d, *J* = 7.8 Hz, 1H), 7.69 (d, *J* = 8.9 Hz, 2H), 7.48 (t, *J* = 7.9 Hz, 1H), 7.26 (d, *J* = 8.5 Hz, 1H), 6.90 (d, *J* = 9.0 Hz, 2H), 3.01 (s, 3H), 2.98 (t, *J* = 5.0 Hz, 4H), 2.84 (t, *J* = 4.9 Hz, 4H).

