

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 8.14 (s, 1H), 8.11 (d, *J* = 8.2 Hz, 1H), 7.68 (s, 1H), 7.62 (ddd, *J* = 11.0, 8.3, 6.3 Hz, 3H), 7.32 (d, *J* = 1.8 Hz, 1H), 7.25 (s, 1H), 7.22 (dd, *J* = 8.2, 1.8 Hz, 1H), 6.05 (s, 2H), 4.26 (t, *J* = 8.4 Hz, 2H), 4.03 (s, 2H), 3.73 (s, 3H), 3.25 (t, *J* = 8.4 Hz, 2H).

