

1H-NMR

¹H NMR (400 MHz, Chloroform-*d*) δ 8.35 (dd, *J* = 8.2, 1.4 Hz, 1H), 8.18 (s, 1H), 7.55 (t, *J* = 7.5 Hz, 1H), 7.34 (d, *J* = 8.3 Hz, 1H), 7.30 (dd, *J* = 8.2, 6.9 Hz, 1H), 2.85 (t, *J* = 7.7 Hz, 2H), 1.78 – 1.74 (m, 2H), 1.42 – 1.25 (m, 8H), 0.88 (t, *J* = 6.7 Hz, 3H).

