

¹H-NMR

¹H NMR (300 MHz, DMSO-*d*₆) δ 10.81 (s, 1H), 7.96 (t, *J* = 5.9 Hz, 1H), 7.49 (d, *J* = 7.5 Hz, 1H), 7.40 – 7.11 (m, 6H), 6.99 (dd, *J* = 16.3, 7.2, 1.3 Hz, 2H), 4.06 (s, 2H), 3.19 (dd, *J* = 8.4, 5.9 Hz, 2H), 2.93 – 2.78 (m, 2H), 1.79 (s, 3H).

