


¹H-NMR

¹H NMR (400 MHz, Deuterium Oxide) δ 5.24 – 5.11 (m, 1H), 5.01 (d, *J* = 4.0 Hz, 1H), 4.08 – 3.91 (m, 2H), 3.74 (dd, *J* = 10.7, 4.0 Hz, 1H), 3.53 (t, *J* = 9.4 Hz, 1H), 3.32 – 3.15 (m, 3H), 2.92 – 2.77 (m, 3H), 2.77 – 2.63 (m, 2H), 2.56 (dd, *J* = 26.6, 10.3 Hz, 1H), 2.46 (s, 3H), 2.38 (s, 2H), 1.94 – 1.84 (m, 1H), 1.76 – 1.55 (m, 3H), 1.44 – 1.31 (m, 1H), 1.14 (s, 4H).



MS

