

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 13.41 (s, 1H), 9.60 (s, 1H), 8.41 (s, 1H), 7.38 (d, *J* = 8.7 Hz, 2H), 6.83 (d, *J* = 8.7 Hz, 2H), 6.69 (d, *J* = 10.1 Hz, 1H), 5.79 (d, *J* = 9.9 Hz, 1H), 5.20 – 5.10 (m, 1H), 3.24 (d, *J* = 7.4 Hz, 2H), 1.78 – 1.72 (m, 3H), 1.65 – 1.61 (m, 3H), 1.44 (s, 6H).

