

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 7.87 (s, 1H), 7.76 (d, *J* = 7.9 Hz, 1H), 7.74 – 7.68 (m, 2H), 7.50 (s, 1H), 7.48 – 7.38 (m, 2H), 7.31 (d, *J* = 7.9 Hz, 1H), 7.19 (ddt, *J* = 8.5, 7.3, 1.1 Hz, 1H), 7.15 – 6.98 (m, 4H), 6.86 (dd, *J* = 16.7, 10.5 Hz, 1H), 6.12 (dd, *J* = 16.7, 2.5 Hz, 1H), 5.68 (dd, *J* = 10.4, 2.5 Hz, 1H), 4.58 (d, *J* = 12.7 Hz, 1H), 4.19 (d, *J* = 13.4 Hz, 1H), 3.21 (t, *J* = 12.7 Hz, 1H), 3.07 (ddd, *J* = 11.8, 8.1, 3.8 Hz, 1H), 2.78 (t, *J* = 12.4 Hz, 1H), 1.94 (d, *J* = 12.9 Hz, 2H), 1.67 (dd, *J* = 22.6, 11.9 Hz, 2H).

