

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 8.85 (t, *J* = 6.1 Hz, 1H), 7.52 (dd, *J* = 9.2, 3.3 Hz, 1H), 7.45 (q, *J* = 8.3 Hz, 4H), 7.34 (ddd, *J* = 9.0, 7.9, 3.3 Hz, 1H), 7.19 (dd, *J* = 9.1, 4.3 Hz, 1H), 6.67 (s, 2H), 5.29 (p, *J* = 7.1 Hz, 1H), 4.56 (d, *J* = 6.1 Hz, 2H), 3.90 (s, 3H), 1.62 (d, *J* = 6.8 Hz, 3H).

