

**<sup>1</sup>H-NMR**

<sup>1</sup>H NMR (400 MHz, DMSO-*d*<sub>6</sub>) δ 11.80 (s, 1H), 8.11 (dd, *J* = 8.1, 1.4 Hz, 1H), 7.91 (d, *J* = 7.4 Hz, 1H), 7.64 (ddd, *J* = 8.5, 6.9, 1.6 Hz, 1H), 7.55 (dd, *J* = 8.5, 1.3 Hz, 1H), 7.31 (ddd, *J* = 8.1, 6.9, 1.2 Hz, 1H), 6.06 (d, *J* = 7.4 Hz, 1H).

