

¹H-NMR

¹H NMR (400 MHz, Chloroform-*d*) δ 6.48 (s, 1H), 4.16 (s, 1H), 2.60 (t, *J* = 6.8 Hz, 2H), 2.10 (s, 3H), 2.09 (s, 3H), 1.86 – 1.71 (m, 2H), 1.55 – 1.46 (m, 3H), 1.44 – 1.32 (m, 4H), 1.32 – 1.20 (m, 10H), 1.19 – 0.99 (m, 7H), 0.87 (s, 3H), 0.86 – 0.82 (m, 9H).

