

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 7.68 (d, *J* = 8.1 Hz, 2H), 7.53 (d, *J* = 7.9 Hz, 2H), 6.91 (d, *J* = 8.6 Hz, 2H), 6.77 (d, *J* = 8.2 Hz, 1H), 6.03 (s, 1H), 3.93 (s, 1H), 3.85 (d, *J* = 10.6 Hz, 1H), 3.74 (s, 1H), 3.24 (d, *J* = 7.1 Hz, 2H), 2.10 (s, 4H), 1.56 (d, *J* = 12.5 Hz, 1H), 1.46 – 1.31 (m, 2H).

