

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 11.65 (s, 1H), 8.10 (d, *J* = 5.5 Hz, 1H), 7.14 (s, 1H), 6.57 (d, *J* = 3.5 Hz, 1H), 4.94 – 4.77 (m, 1H), 4.20 – 4.01 (m, 2H), 4.00 – 3.88 (m, 1H), 3.88 – 3.60 (m, 2H), 3.44 – 3.39 (m, 1H), 3.27 (d, *J* = 4.0 Hz, 3H), 2.46 – 2.32 (m, 1H), 1.90 – 1.66 (m, 1H), 1.65 – 1.51 (m, 1H), 1.01 (d, *J* = 7.2 Hz, 3H).

