

¹H-NMR

¹H NMR (400 MHz, DMSO-*d*₆) δ 9.58 (s, 1H), 9.06 (s, 1H), 8.40 (s, 1H), 8.38 – 8.33 (m, 1H), 8.09 (d, *J* = 5.3 Hz, 2H), 7.40 – 7.29 (m, 2H), 7.20 – 7.10 (m, 2H), 4.38 (d, *J* = 5.7 Hz, 2H), 3.77 – 3.71 (m, 1H), 3.58 – 3.50 (m, 3H), 3.20 – 3.13 (m, 1H), 2.82 (d, *J* = 13.7 Hz, 1H), 1.76 (q, *J* = 7.4, 6.8 Hz, 2H), 1.65 – 1.06 (m, 14H).

