

产品说明书

26 版

基本信息

产品编号	产品名称	质量标准	粒径 (μm)	电荷	pH	比表面积 (m ² /g)
C80129	AM 胶态二氧化硅	30 wt. % suspension in H2O	12	负	9	220
C80130	AS-30 胶态二氧化硅	30 wt. % suspension in H2O	12	负	9.4	228
C80131	AS-40 胶态二氧化硅	40 wt. % suspension in H2O	22	负	9.4	228
C80132	CL 胶态二氧化硅	30 wt. % suspension in H2O	12	正	5	230
C80133	CL-P 胶态二氧化硅	41 wt. % suspension in H2O	22	正	4	140
C80134	CL-X 胶态二氧化硅	46 wt. % suspension in H2O	22	负	9	130
C80135	DF 胶态二氧化硅	30 wt. % suspension in H2O	12	负	9.8	228
C80136	HS-30 胶态二氧化硅	30 wt. % suspension in H2O	12	负	9.8	228
C80137	HS-40 胶态二氧化硅	40 wt. % suspension in H2O	12	负	9.5	228
C80138	LS 胶态二氧化硅	30 wt. % suspension in H2O	12	负	8.2	228
C80139	PG-E 胶态二氧化硅	30 wt. % suspension in H2O	22	负	9	155
C80140	SK 胶态二氧化硅	25 wt. % suspension in H2O	12	负	4.2	
C80141	SM 胶态二氧化硅	30 wt. % suspension in H2O	7	负	10	360
C80142	SM-30 胶态二氧化硅	30 wt. % suspension in H2O	7	负	10	360
C80143	TM-40 胶态二氧化硅	40 wt. % suspension in H2O	22	负	9	130
C80144	TM-50 胶态二氧化硅	50 wt. % suspension in H2O	22	负	9	130
C80145	TMA 胶态二氧化硅	34 wt. % suspension in H2O	22	负	7	130

储存: 室温

