



## Specification Sheet

|                   |                                                                              |
|-------------------|------------------------------------------------------------------------------|
| Product Name:     | T0080                                                                        |
| Product Number:   | T14738                                                                       |
| CAS:              | 2785323-68-2                                                                 |
| Formula:          | C <sub>24</sub> H <sub>22</sub> F <sub>3</sub> N <sub>3</sub> O <sub>3</sub> |
| Molecular Weight: | 457.44                                                                       |
| Storage at:       | -20°C                                                                        |

| Test                                     | Specification         |
|------------------------------------------|-----------------------|
| Appearance                               | White Powder          |
| Purity(HPLC)                             | ≥98%                  |
| <sup>1</sup> H-NMR(DMSO-d <sub>6</sub> ) | Conforms to Structure |

